MODEL

PROFESSIONAL 2.0 USER MANUAL

ENG

CZ

SK

LANGUAGES:

ER 01

alcoforce.com

BREATHALYZER
ELECTROCHEMICAL

Dear Customer, congratulations on your choice! We invite you to read the user's manual, so that the purchased device serves you in best possible way.

SCHEME & SPECS

SENSOR
platinum electrochemical

PRECISION
 $\pm 0.06\%$ at 0.5%
in 25°C

TESTING RANGE
 $0,000\%$ - $3,000\%$

RE-CALIBRATION
every 12 months
or 600 tests

SIZE AND WEIGHT
 $46 \times 115 \times 27\text{mm}$,
 $77,5\text{g}$ (with batteries)

POWER SUPPLY
1 x AAA 1,5V
(LR03) alkaline

MOUTHPIECE
required

OTHER
The display turns red when
the result of the measurement
exceeds the threshold
of intoxication

measuring hole

display

function
button

mouthpiece

air outlet


battery
cover

PREPARING THE DEVICE FOR WORK


1

Remove the back battery flap by pressing the triangle on the flap and slide it downward.

2

Place **1.5V AAA LR03 alkaline batteries** in the power compartment according to polarity (plus up, minus down). The strap must be under the bottom of the battery (it is used to pull the battery out). After placing the batteries in the power compartment, press the battery flap against the device guides and push it up until you hear a click.

* do not use rechargeable batteries!

STORAGE AND GENERAL REMARK

30 min

Do not eat, drink or smoke for 30 minutes before taking the measurement.

0-40 °C

Store in a dry place, avoid direct exposure to sunlight and temperatures outside the range of 0 to 40 °C. Do not store in the car!

Do not take a measurement directly after consuming alcohol. Blowing into the breathalyzer immediately after consuming alcohol can damage the sensor and cause incorrect measurements.

Do not store the breathalyzer near volatile substances including alcohol. Do not blow liquids or fumes into the sensor or immerse the device in liquids. Do not clean the device with detergents.

Service your device regularly or whenever it displays a service message. Within 12 months from the date of purchase, service is completely free* - take advantage of this. Avoiding service will negatively affect the life and accuracy of the breathalyzer!

*NOTE. Applies to purchase by the first user.

Do not disinfect or clean the device with alcohols or other disinfectants. Do not touch the breathalyzer for at least 30 minutes after hand disinfection.

ENG

Rinse your mouth after drinking the alcohol and wait 30 minutes. This amount of time is crucial for the alcohol to absorb into your bloodstream and to reduce the alcohol vapours in your mouth and saliva. **For the correct results, do not eat, drink or smoke cigarettes for min. 30 minutes before taking the measurement.** Do not take measurements in strong winds, in air-conditioned rooms, or in polluted air. Do not use immediately after hand disinfection or in the rooms with present disinfectants.

Many factors may interfere with the measurement results (heavy environmental pollution, smoking, various gas vapors, wind, air conditioning, alcohol vapors coming from disinfectants), hence the results presented by the breathalyzer are only estimates. The manufacturer, importer and seller do not bear any responsibility for the correctness of the results and potential damage caused by improper use of the breathalyzer. The device cannot be used as an evidence breathalyzer.

Follow the information on the breathalyzer's display and in the instructions when taking a measurement.

30 minutes

End smoking,
eating and drinking

Rinse your mouth


Switch on
the device

Take 3
measurements at
1 minute intervals.
Choose the highest
score.

HOURLY ALCOHOL ELIMINATION

The alcohol level in blood and exhaled air changes over time. At first, it increases during the absorption of alcohol into the blood, then decreases during its elimination (burning) in the body. A breathalyzer can help you estimate which phase you are in: alcohol concentration increase or elimination. To obtain a reliable result, take measurements when you no longer intend to consume alcohol. We suggest to make min. 3 measurements at equal 30 minute intervals and to record the individual results. Comparison of these results will give an answer to the question in which phase you are - using 3 decreasing results obtained in equal intervals it is possible to estimate when the body will get rid of alcohol.

3 decreasing measurements taken at intervals of 30 minutes.

CARRYING OUT THE MEASUREMENT

1

Insert the mouthpiece into the measuring hole (*1).

Turn on the device by pressing the function button on the main panel (*2).

2

1

2

3

After starting, you will hear a single beep. The display will show four eights (*1), and then the breathalyzer will display the number of measurements already taken (*2) and start preparing for the test by counting down from 10 to 0 (*3).

3

The device is ready to take an air sample. Take a deep breath and blow into the mouthpiece for approx. 6 sec. The device will emit an intermittent beep while the sample is being taken.

4

The beep ended. The breathalyzer has analyzed the sample and displayed the result. If the result is above 0.20‰, the display will light up red. After 10 s. the result will disappear and the breathalyzer will automatically turn off.

sample results

ENC

MESSAGES

Provided air sample is too small. Perform the measurement again providing a strong and steady exhalation.

Low voltage. Replace the batteries.

The measurement result exceeds the measurement range of the device.

TROUBLESHOOTING: READ BEFORE YOU RETURN FOR REPAIR

If you have the impression that your breathalyzer is not working properly, we recommend to check the following:

1. The use of old or improper batteries may affect the accuracy of the results. Use only high-quality alkaline batteries. **Do not use rechargeable batteries.**
2. Decalibration of the breathalyzer. Each breathalyzer should be calibrated in accordance with its technical specification. For 12 months from the purchase of the breathalyzer, we will perform the calibration service for **FREE!!!** For more information, please visit alcoforce.com.
3. Each person has a different ability to metabolize alcohol, and this ability can change periodically depending on a number of factors. So if your **PROFESSIONAL 2.0** Breathalyzer is regularly calibrated according to the service interval and you use it according to the instructions and the following three results are close to each other, it is most likely working properly. You can read more about the factors that affect breathalyzer results at: alcoforce.com.

THE SELLER DOES NOT ACCEPT SERVICE REQUESTS


service requests are considered by the central service
AlcoForce / Work and Road Safety Corporation

Calibration is a concept inherent in EVERY Breathalyzer.

This means that every 12 months or so the breathalyzer should be reconfigured to correctly interpret the indications of the wearable sensor. The need for calibration is evidenced by the messages displayed on the device and results that deviate significantly from reality - for example, 0.5‰ despite the absence of alcohol consumption.


The necessity of calibration is not a disadvantage of the breathalyzer and is not subject to complaint. Please do not expect that a breathalyzer that is not calibrated regularly will show the correct results. Calibrations after the free service period are chargeable and can be made at the premises of ALCOFORCE / WORK AND ROAD SAFETY CORPORATION (also by mail order) or at another authorized breathalyzer service point. **Additional information can be found in the Warranty Card.**

Calibration in this model should be performed every **600 measurement** or **12 months** (regardless of the number of measurements performed). Regular calibration is necessary to maintain proper operation of the device and extends the life of the sensor.

12 months

free calibration

Note - calibration of this breathalyzer within **12 months** from the date of purchase is **FREE OF CHARGE** at the premises of AlcoForce and Work and Road Safety Corporation of other authorized service points. **More information can be found in the Warranty Card.**

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