

EV109

GB	Digital personal scale with BMI index
CZ	Digitální osobní váha s BMI indexem
SK	Digitálna osobná váha s BMI indexom
PL	Cyfrowa waga łazienkowa z pomiarem BMI
HU	Digitális személymérleg (BMI)
SI	Digitalna osebna tehtnica z indeksom telesne mase
RS HR BA ME	Digitalna lična vaga sa BMI indeksom
DE	Digitale Personenwaage mit BMI-Index
UA	Цифрові особисті ваги з BMI показником
RO	Cântar digital de persoane cu indice BMI
LT	Skaitmeninės svarstyklės at BMI indeksu
LV	Digitālie personiskie svāri su BMI indeksu
EE	Digitaalne kaal BMI indeksiga
BG	Цифрова везна за баня с индекс на BMI

1

2

GB | Digital personal scale with BMI index

Read this manual and safety instructions carefully before first use of the product!

Technical specifications:

Toughened safety glass: 6 mm
 LCD display: 91 × 43 mm
 Maximum weight: 180 kg
 Minimum load: 5 kg
 Measurement increments: 0.1 kg / 0.2 lb
 Unit of measure: kg/lb/st
 Memory: max. 13 users
 Height setting: 100–240 cm

Age setting: 10 to 99 years
 Unit, increment of fat content measurement: 0.1 %
 Unit, increment of water content measurement: 0.1 %
 Unit, increment of muscle mass measurement: 0.1 %
 Automatic switch-off after 10 seconds
 Indicator for low battery/overload
 Power supply: 3× 1.5 V AAA (included)

Basic rules for use:

1. Digital bathroom scales EV109 are designed for household use only.
2. Have all repairs done by a professional in a repair shop or the store where you bought the device!
3. Do not use the scales in the vicinity of devices which create electromagnetic fields. Place the scales at a sufficient distance from such products to prevent mutual interference.
4. The result is only an approximate value. The product is not a medical tool. Please always consult your doctor regarding medication or diet to achieve optimal weight.
5. Do not insert any objects into the openings on the device.
6. The scales are intended for weighing persons older than 10 and younger than 99 years of age.
7. Weighing house pets is not allowed.
8. Protect the scales from falls or impacts. Do not jump on the scales.
9. You or the device may slip when on a wet surface. Similarly, never stand on the edge of the scales when on a smooth surface, and keep balance while standing on the scales.
10. Due to electronic impedance, never use the device in the following circumstances:
 - a) if you have an implanted pacemaker or another medical electronic implant
 - b) if you suffer from heart problems or disease
 - c) during pregnancy
 - d) if you are treated with dialysis
 - e) if you have a fever
 - f) if you suffer from osteoporosis
 - g) if you have an oedema or swelling
 - h) if you have been exercising longer than 5 hours
 - i) if your heart rate is lower than 60 beats per minute
11. Only use the scales in accordance with the instructions provided in this manual.
12. The producer is not liable for damage caused by improper use of the device.
13. This appliance is not intended for use by persons (including children) whose physical, sensory or mental disability or lack of experience and expertise prevents safe use, unless they are supervised or instructed in the use of the appliance by a person responsible for their safety. Children must always be supervised and must never play with the appliance.

Button description:

(see Fig. 1)

- 1 – ▲ button – increases value by 1 unit of measure
- 2 – ■ button – manually turns on the scales, input information settings
- 3 – ▼ button – reduces value by 1 unit of measure

Back

(see Fig. 2)

- 1 – button for setting mass units
- 2 – battery compartment

Overview of Icons

Icon	Meaning
kg st:lb	Body mass unit
P01	User position in memory
	Male
	Female
cm	Body height
Age	Age
FAT	Body fat
TBW	Water content in body
	Muscle mass
	Bone mass
%	Values in percent – body fat / water / muscle mass / bone mass
Underfat Healthy Overfat Obese	BMI rating – underweight / normal / overweight / obese
KCAL	Daily kilocalorie requirement
BMI	BMI value

Inserting/Replacing Batteries

1. Open the battery compartment cover on the bottom side of the scales.
2. Remove the depleted batteries.
3. Insert new 3x 1.5 V AAA batteries. Make sure you observe the correct polarity.
If batteries are protected by a safety film, remove it first before putting the batteries into the scales.
Only use alkaline batteries of the same type; do not use old and new batteries at the same time; do not use rechargeable batteries.
4. Close the cover.

Basic Weight Measurement (without BMI) / Changing Unit of Weight

1. Place the scales on a flat, hard and stable surface.
2. Carefully step onto the scales - the scales turn on automatically.
If you wish to change the unit of mass (kg/lb/st), repeatedly press the mass unit change button on the underside of the scales.
3. Stand still on the scales and do not move – after a moment, the measured weight will be displayed.
4. After 10 seconds, the scales turn off automatically.

MEASURING BODY FAT, WATER, MUSCLE MASS, BONES

INFORMATION FOR MEASURING BODY FAT

WARNING: Persons using a pacemaker or another implant with electronic components are not allowed to use these scales – proper functioning of the implant might be impaired.

NOTE: Correct measurements can only be achieved when standing evenly on the scales, barefoot and with dry feet! The basic principle of the device is measuring the electrical impedance within the human body. To obtain data, the scales send a weak, unnoticeable electronic signal through the body.

The signal is completely safe! The method is called „Bioelectrical Impedance Analysis“ (BIA), a factor based on the mutual relation between the body and weight and the associated percentage of water in the body, as well as other biological data (age, sex, body mass). There are corpulent people with a balanced ratio between fat and body mass and lean people with very high body fat value. Body fat is a factor which can significantly impact for instance the cardiovascular system. Therefore, checking both values (fat and weight) is a useful tool for maintaining good fitness and a healthy body.

BODY FAT

Body fat is an important component of the body – from protecting the joints, to being a source of vitamins, a body temperature regulator, and overall a vital role player in the human body. It is advised to not reduce its amount rapidly, but to ensure it is at a healthy level, keeping a good balance between muscle mass and fat.

IMPORTANCE OF WATER FOR OUR HEALTH

The human body is approximately 55–60 % water, with ratios varying based on age and sex.

Water has several roles in the human body:

- constitutes our cell blocks; meaning that every cell, such as skin cells, brain cells, muscle cells etc. operate correctly only if they contain sufficient water,
- it is an important substance released in our body,
- it transports nutrients and metabolites in the body.

You can compare the reference values for your fitness profile with the following table. It shows the dependence of body fat on age and sex. Focus more on keeping a balance of water in the body. Similarly to fat levels, water levels fluctuate during the day. The obtain relevant results, make measurements always at the same time of day. The obtained measurements for body fat and water should be taken rather as informative. They are not recommended for medical purposes.

Always consult these measurements with your doctor before making radical changes to your regimen.

Measuring Weight with Display of Body Fat, Water, Muscle Mass, Bones, BMI

Before weighing, enter basic information about each user.

Can save data for up to 13 users.

Take off your shoes and socks before weighing, otherwise the information displayed will be incorrect.

Saving User Information

Press the ■ button.

P01 will begin flashing on the screen, indicating user no. 1.

Again press the ■ button to save the user's information under number P01 or press the Λ/V button to change the user number.

Input entry information using buttons Λ/V in the following order:

User number – sex – age – height.

Confirm each value by pressing ■.

You can later change the saved information by pressing the ■ button.

After inputting entry information, wait until the display shows 0.0 kg.

Step onto the scales and stand still; do not move – after a moment, the display will show:

Current weight → body fat → water → muscle → bones → daily kilocalorie requirement → BMI.

This information will be displayed 2x in a row.

Then the scales will turn off.

BMI

BMI	Underweight	Normal	Overweight	Obese
	<18.5	18.5–25	25–30	>30

Amount of Fat in the Body

BODY FAT – MALE

AGE	Underweight	Normal	Overweight	Obesity
<20	<19	19–29	29–35	>35
20–30	<20	20–30	30–36	>36
30–40	<21	21–31	31–37	>37
40–50	<22	22–32	32–38	>38
50–60	<23	23–33	33–39	>39
>60	<24	24–34	34–40	>40

BODY FAT – FEMALE

AGE	Underweight	Normal	Overweight	Obesity
<20	<23	23–33	33–39	>39
20–30	<24	24–34	34–40	>40
30–40	<25	25–35	35–41	>41
40–50	<26	26–36	36–42	>42
50–60	<27	27–37	37–43	>43
>60	<28	28–38	38–44	>44

The scales may show the following error messages:

Lo – low batteries, replace.

Err – scales overloaded, step off the scales. The scales could become damaged.

ErrL – fat level too low.

Errh – fat level too high.

C – error during measurement, reweigh.

Do not dispose of electrical appliances as unsorted municipal waste; use sorted waste collection points. Contact local authorities for up-to-date information about collection points. If electrical appliances are deposited in waste landfills, hazardous substances may leak into the groundwater, enter the food chain and harm your health.

Emos spol s.r.o. declares that the EV109 is in compliance with the essential requirements and other applicable provisions. The devices can be freely operated in the EU. The Declaration of Conformity can be found at <http://www.emos.eu/download>.

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